





## an-Lens UPGRADE

### Take Advantage of the MD8 Design



The **an-lens MD6** and **MD8** have the:

- **same optic size**
- **same haptic design**
- **same optical design**

The **MD8** offers the added advantage of:

- **Additional 360° barrier to delay onset of PCO**

#### Use of the MD8

- To take advantage of the posterior barrier, simply load the lens inside the cartridge using the clearly marked notches on the haptics
- If the lens is implanted the other way, it remains identical to the MD6

We offer the advanced MD8 for the **same price** as the MD6!